

DK Ramper 3
GB/US	... Ramps 6
DE Rampen 9
SE Ramper 12
FR Rampes transportables 15
IT Rampe 18
NO Rampe Skinne/ramper 21
FI Luiskat. 24
U.S.A. and countries outside the EU	27


Version 5.00


I Fast rampe og bred fast rampe
Plain ramp and wide plain ramp
Feste Rampe und breite feste Rampe
Fast ramp och bred fast ramp
Rampe étroite fixe et rampe large fixe
Rampa fissa e rampa larga fissa
Skinne fast og bred fast rampe
Kiinteä tasoluiska ja leveä tasoluiska


II Længdeforskydelig rampe
Telescopic ramp
Längenverschiebbare Rampe
Längdjusterbar ramp
Rampe télescopique
Rampa telescopica
Skinne teleskop
Teleskoopluiska

III Sammenklappelig og bred sammenklappelig rampe
Folding ramp and wide folding ramps
Zusammenklappbare Rampe und breite zusammenklappbare Rampe
Ihopfällbar ramp och bred ihopfällbar ramp
Rampe étroite pliante et rampe large pliante
Rampa pieghevole e rampa larga pieghevole
Skinne foldbar og rampe bred foldbar
Taittoluiska ja leveä taittoluiska


IV EasyFold-Pro

V EasyFold-Pro3

VI Længdeforskydelig kuffertrampe
EasyFold Telescopic ramp
Längenverschiebbare Kofferrampe
Längdjusterbar och vikbar ramp
Rampe portefeuille télescopique
Rampa a valigetta telescopica
Rampe koffert foldbar, teleskop
Salkkuluiska teleskoopilla

Varenumre

30000 - xxx
 30010 - xxx
 30015 - xxx
 30020 - xxx
 30030 - xxx
 30040 - xxx
 30070 - xxx
 30080 - xxx
 30085 - xxx

Contents

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1.00	Purpose and Application	6
1.01	Manufacturer	6
1.02	Intended Purpose	6
1.03	Conditions for Use	6
1.04	Important/Warnings	7
1.05	Unpacking and Preparation of the Ramp	7
2.00	Description of Functions	7
2.01	Pictograms/Labels	7
3.00	Transportation	7
4.00	Maintenance and Storage	7
4.01	Cleaning the ramps	7
4.02	How to Prevent/Avoid Corrosion	7
4.03	What Maintenance Is Required By the Owner?	7
5.00	Service and Useful Life	8
5.01	Useful Life and Servicing of the Ramp	8
5.02	Inspection of Stepless ramps	8
6.00	Technical Specifications	8
7.00	EC-Declaration of Conformity	8

1.00 Purpose and Application

1.01 Manufacturer

V. Guldmann A/S
 Graham Bells Vej 21-23A
 DK - 8200 Aarhus N
 Tel. +45 8741 3100

1.02 Intended Purpose

Stepless Ramps are intended for helping wheelchair users overcome height differences in the physical surroundings.


Stepless ramps are a carefully designed and fully developed series of transportable wheelchair ramps which can be used where wheelchair users need to overcome height differences of up to 60 cm.

Stepless ramps are available in four basic versions:
 Wide ramps, available in both fixed and folding models.
 Narrow ramps, available in fixed, telescopic and folding models.
 Wide ramps available as telescopic suitcase ramps and EasyFold-Pro.
 EasFold-Pro3 that folds in three parts.

1.03 Conditions for Use

Placement of EasyFold-Pro and EasyFold-Pro3

Place the ramp correctly before usage, open the ramp after placement. See pictures for correct placement.

Plain ramp, wide plain ramp, EasyFold-Pro and EasyFold-Pro3

The ramp should be placed on a firm underlay. When setting up narrow ramps, the two parts of the ramp should be placed parallel to one another at a distance corresponding to the distance between the wheels of the wheelchair. The ramp should be positioned such that the entire bent contact surface of the ramp is supported. When setting up EasyFold-Pro, the ramp should be positioned such that the bent contact surface is supported. Then unfold the upper ramp. EasyFold-Pro3 folds in three parts, so it is easier to store away. To unfold the ramp, pull the lock beside the handle and unfold the ramp.

Telescopic ramp

The ramp should be placed on a firm underlay. When setting up narrow ramps, the two parts of the ramp should be placed parallel to one another at a distance corresponding to the distance between the wheels of the wheelchair. The ramp should be positioned such that the entire bent contact surface of the ramp is supported.

To extend the ramp

Hold the ramp by the handle in your right hand, and pull the ramp out to its full length by gripping the contact surface of the inner part of the ramp. Position as required. If the catch device should be moving too tightly, it can be lubricated with oil or another lubricant.


Telescopic suitcase ramp

The ramp should be placed on a firm underlay.

The ramp should be positioned such that the entire bent contact surface of the ramp is supported.

The ramp should be carried by the handles and the contact surface positioned on the upper level. Then unfold the upper ramp. The upper ramp closest to the body is extended to the desired length and placed on the upper level. Afterwards the other half of the upper ramp is extended. Always check that the ramp is placed in a stable and safe way.

To prevent the suitcase ramps from sliding on underlays with low friction, the safety pin should always be mounted before using the ramps.

Important!

Do not start up electric wheelchairs which have stopped on the ramp, as the various parts of the ramp may extend unintentionally by so doing. In the event of a wheelchair stopping in the middle of the ramp, reverse down and start from the bottom again.

Folding ramp, wide folding ramps

When setting up narrow ramps, the two parts of the ramp should be placed parallel to one another at a distance corresponding to the distance between the wheels of the wheelchair.

The ramp should be positioned such that the entire bent contact surface of the ramp is supported.

The ramp should be placed on a firm underlay.

Important!

Unfold the ramp by gripping both handles.

To avoid anything getting caught up in the hinge, the ramp should be held with the hinge opening away from the body whilst it is being opened and closed.

1.04 Important/Warnings

- Read all the operating instructions before using the ramp.
- Do not exceed the load marked on the ramp.
- The ramp should only be used for driving wheelchairs up and down.
- Wheelchairs should be driven up and down at low speed
 - without a run-up.
- The ramp should be positioned such that the entire contact surface is supported.
- The ramp should be used on a firm underlay.
- The ramp should not be used for gradients greater than 1:5
 - for electric wheelchairs or wheelchairs pushed by a helper.
- The ramp should not be used for gradients greater than 1:7
 - independent user of manual wheelchair.

- Watch out not to squeeze fingers when closing the ramp.
- The handle on EasyFold-Pro can not be changed if it breaks off, it will only break off if it is not placed back in place after use.
- Do not kick or push the ramp in place
- If damage is suspected, the ramp must not be used before the dealer has authorized it.
- Any serious incident that occurred in relation to this device should be reported to the manufacturer and the local competent authority.

1.05 Unpacking and Preparation of the Ramp

Visual Check of the Ramp

If the packaging is damaged on receipt, the ramp should be inspected carefully for visible deficiencies, defects or faults.

If damage is suspected, the ramp must not be used before qualified service personnel or the Guldmann Service Team have authorised it.

2.00 Description of Functions

2.01 Pictograms/Labels

CE-marking

- Serial number label
- Stepless label
- Label stating max. load

3.00 Transportation

The ramp should be carried by the handle(s).

4.00 Maintenance and Storage

4.01 Cleaning the ramps

Clean the ramps with a damp cloth, sponge or soft brush applied with a normal cleaning agent.

The ramps may also be rinsed under running water from garden hoses or the like.

Chemicals and autoclaving may not be used to clean the ramp.

4.02 How to Prevent/Avoid Corrosion

The ramp is made from corrosion-resistant materials, and does not require further corrosion prevention.

4.03 What Maintenance Is Required By the Owner?

The hinges on the Stepless folding ramps should be lubricated with TRI-FLOW or a similar lubricating agent every six months, or as required.

The sliding rails of the Stepless telescopic ramps should be dried off with a dry cloth when needed.

If the ramps are slippery (covered with ice), use salt on the ramp - NEVER gravel.

The telescope system must not be taken apart by the user

If the telescope system of the telescopic ramp is taken apart, it should be taken apart and put together again by qualified personnel.

Tightening moments: Hinge = 10 Nm. Other = 5 Nm.

Safety: All screws and nuts should be fastened with Locktite no. 270.

5.00 Service and Useful Life

5.01 Useful Life and Servicing of the Ramp

Given the maintenance prescribed above, Stepless ramps will have a useful life of 15 years.

Spare part lists and drawings can be ordered from the manufacturer or the distributor.

5.02 Inspection of Stepless ramps

Parts which are worn or defective should be replaced by new spare parts from V. Guldmann A/S.

1. Visual check of the ramp

- Check the ramp for wear and tear.
- Check the ramp for any faults.
- Check that there is no other form of damage to the ramp.

2. Testing of the ramp, as for normal use

- Check all the ramp's functions, loaded and unloaded.

3. Control of the mechanical state of the ramp

- Clean dirt and other impurities from the ramp.
- Inspect and assess vital parts of the ramp.
- Replace any parts of the ramp which are defective or worn.

4. Go through point 2 again to recheck

5. Did any new problems arise/were new problems noted in point 4?

- If new problems arose, go back to point 3.
- If no new problems were noted, conclude the inspection.

6.00 Technical Specifications

Plain, telescopic and folding ramps

Lengths 500 - 3100 mm
Widths
outer, max 260 and 279 mm
inner 182 mm
Max. load under 2100/over 2100 mm each pair 350/200 kg

Wide plain and wide folding ramps

Lengths 1000 - 3000 mm
Widths
outer, max 880 mm
inner 770 mm
Max. load under 2100/over 2100 mm 325/200 kg

EasyFold-Pro

Lengths 900 - 2100 mm
Widths
inner 760 mm
folded, max. 398 mm
unfolded, max. 776 mm
Max. load 350 kg

EasyFold-Pro3

Lengths 1635 mm
Widths
inner 760 mm
folded max. 273 mm
unfolded max. 786 mm
Max. load 350 kg

Telescopic suitcase ramp

Lengths 1830 - 3050 mm
inner 760 mm
folded, max. 405 mm
unfolded max. 820 mm
Max. load under 2440/over 2440 mm 325/200 kg

All types

Height of edge 46 mm
Angle of contact max. 12 degrees

Weight, min. and max. per ramp, depending on length

Plain ramp 1.5 - 6.0 kg
Telescopic ramp 3.5 - 10.1 kg
Folding ramp 4.0 - 10.1 kg
Wide plain ramp 9.7 - 28.6 kg
Wide folding ramp 17.5 - 31.7 kg
Telescopic suitcase ramp 16.4 - 26.4 kg
EasyFold-Pro 7,8 - 15,8 kg
EasyFold-Pro3 13,4 kg

Materials

Covering and handle Impact-proof and non-skid plastic
Ramp slopes Anodised aluminium, with skid resistant surface

7.00 EC-Declaration of Conformity

The product is manufactured in compliance with regulation (EU) 2017/745 of the European parliament and of the Council of 5 April 2017, as medical device class 1.

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U.S.A. and countries outside the EU

A. Users guide

Before using the product, read the entire operation manual including warranty..

B. WARRANTY

Guldmann warrants its equipment is free from material defects under normal use, and will perform substantially in accordance with the specifications set forth in documentation provided with the equipment.

This express warranty shall be in effect for one year from the date of original purchase and installation (the "Warranty Period"). If a valid claim is made during the Warranty Period for malfunction or equipment defect, Guldmann will repair or replace the equipment at no additional cost to you. Guldmann retains sole discretion as to whether the equipment will be repaired or replaced.

This warranty shall be null and void if the equipment is operated and maintained in any manner inconsistent with its intended use or the instructions provided with the product. Further, in order for the warranty to remain in effect for the full Warranty Period, all service to the equipment must be provided by a Guldmann designated technician. Any parts or components repaired or replaced by a Guldmann designated technician will be guaranteed for the remainder of the Warranty Period.

The warranty does not cover any part of the equipment which has been subject to damage or abuse by the user or others. The warranty does not cover any part of the equipment which has been altered or changed in any way by the user or others. Guldmann does not warrant that the lifting device functions will meet your requirements, be uninterrupted or error free.

The warranty set forth is in lieu of all other express and implied warranties, whether oral, written or implied, and the remedies set forth above are your sole and exclusive remedies. Only an authorized officer of Guldmann may make modifications to this warranty, or additional warranties binding on Guldmann. Accordingly, additional statements such as advertising or presentations, whether oral or written, do not constitute warranties by Guldmann.

Service or Repair

Contact Guldmann Repair for an authorization to return any defective item during the Warranty Period. You will be provided with a return authorization number and address for returning the item for warranty service or replacement. Do not return items to Guldmann under warranty without receiving a Return Authorization Number.

If mailing the item, pack it carefully in a sturdy carton to prevent damage. Include your Return Authorization Number, a brief description of the problem and your return address and phone number. Guldmann does not assume the risk of loss or damage while in transit, so it is recommended you insure the package.

V. Guldmann A/S
Head Office:
Tel. +45 8741 3100
info@guldmann.com
www.stepless.com

Guldmann GmbH
Tel. +49 611 974 530
Fax +49 611 300 018
info@guldmann.de
www.stepless.de

Guldmann Sverige AB
Tel. +46 0322 55290
info@guldmann.se
www.stepless.se

Guldmann Sarl
Tél. +33 145 54 78 36
Fax +33 145 57 10 52
france@guldmann.com
www.stepless.fr

Guldmann Srl
N° Verde 800 781 604
Tel. +39 0521 660132
Fax +39 0521 691399
Cell. +39 340 9087107
italia@guldmann.com
www.stepless.it

Guldmann Inc.
Tel. 800 664 8834
Tel. 813 880 0619
Fax 813 880 9558
info@guldmann.net
www.stepless.net