



Set Date & Time – One2Clean

With unit turned on			Cycle pause rest: 684:44min
Press			
Press		X5	19:33:40 2014-05-21 we
Press			1□:33:59 2014-05-21
Press		To change hours to correct time	
Press		To move to minutes	09:3□:59 2014-05-21
Press		To change minutes to correct time	
Press		To move to date	09:11:59 201□-05-21
Press		To change the year	
Press		To move to month	09:11:59 2016:0□:21
Press		To change the month	
Press		To move to day	09:11:59 2016-01-2□
Press		To change the day	
Press		To store new settings	
Press		To go back to auto mode	