



USER COMFORT

The Ropox shower couch has soft padded waterproof removable slats for supreme comfort and easy cleaning and maintenance.



CONVENIENT

The Ropox changing couch folds away against the wall to create space when not in use.



SIMPLE OPERATION

An optional adjustable angle backrest section provides comfort and support for users who are unable to lie prone.

ROPOX SHOWER COUCH

A stylish changing solution

With its clean lines and slimline design, the Ropox shower couch with electric operation offers versatility and value for money in a stylish package. The cushioned mattress section gives users the ultimate in comfort, and optional see through side and end safety rails create a safe environment.

A water collection tray and drainage hose ensure safe and efficient removal of waste water, and the easy action folding mechanism enables users to fold the couch away after use, creating space in the room. With a travel length of 70 cm the Ropox shower couch provides safe hygiene for users of all sizes.

It is available in three different lengths; 146, 178 and 210 cm enabling it to be fitted in most bathrooms.

3 LENGTHS

ADJUSTABLE

EASY TO CLEAN




ROPOX SHOWER COUCH



FOLDED DEPTH
314mm



OVERALL DEPTH
942mm



SURFACE DEPTH
760mm



MIN HEIGHT
300mm



MAX HEIGHT
1000mm



HEIGHT RANGE
700mm



LENGTH
1410mm, 1910mm,
2100mm



SWL
220kg



WARRANTY
1 year



SERVICES

Installation, commissioning,
maintenance, servicing,
training

Call 0800 298 6000 today to arrange an assessment for the Ropox Couch

OPTIONS & ACCESSORIES

LENGTH WITH/WITHOUT BACK SUPPORT 146cm, 178cm, 210cm

SAFETY SIDE RAIL for user comfort & reduced entrapment risk

SAFETY END RAIL for additional security

BACKREST for additional comfort & support, to raise users to a semi seated position

www.multicaremedical.co.uk

0800 298 6000

sales@multicaremedical.co.uk



Multicare Medical Ltd, Progress House, Jacknell Road, Dodwells Bridge Industrial Estate,
Hinckley, Leicestershire, LE10 3BS

V01